
FITTING YOUR MASK

F&P Brevida™

NASAL PILLOWS MASK

Use and Care Guide

F&P Brevida™

**Thank you for choosing
your F&P Brevida.**


This nasal pillows mask has been designed to make it easier for you to fit CPAP therapy into your life. This is your mask Use and Care Guide. Please read this entire guide before you first use your mask. Following the steps in this guide, and referring to the images and diagrams will assist you with assembly, fitting and daily use of your F&P Brevida nasal pillows mask.

REF 185048407 REV D 2019-02 © 2019 Fisher & Paykel Healthcare Limited

Fisher & Paykel
HEALTHCARE

A

Fisher & Paykel
HEALTHCARE

Rx only

CE 0123

LANGUAGES

PAGE	LANGUAGE
1	English
13	Français (French)
25	Español (Spanish)
37	Bahasa Melayu (Malaysian)
49	ไทย (Thai)
61	简体中文版 (Simplified Chinese)
73	繁體中文版 (Traditional Chinese)
85	한국어 (Korean)
97	日本語 (Japanese)

Manufacturer
Fisher & Paykel Healthcare Ltd,
15 Maurice Paykel Place,
East Tamaki, Auckland 2013
New Zealand
Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Importer/ Distributor
Australia (Sponsor)
Fisher & Paykel Healthcare Pty Ltd,
19-31 King Street, Nunawading,
Melbourne, Victoria 3131.
Tel: +61 3 9871 4900
Fax: +61 3 9871 4998

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USA/Canada
Tel: 1800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001

INTENDED USE

The F&P Brevida Nasal Pillows Mask is intended to be used by individuals who have been diagnosed by a physician as requiring CPAP or Bi-Level therapy. The F&P Brevida Nasal Pillows Mask is intended for single patient adult (≥ 66 lbs (30 kg)) use in the home and for multiple patient adult use in the hospital or other clinical setting where proper disinfection of the device can occur between patient uses.

OPERATING INSTRUCTIONS

- The operating pressure range of your mask is 4 to 25 cmH₂O.
- The operating temperature range of your mask is 5 to 40 °C (41 to 104 °F).
- Before using your mask each time, you should:
 - i. Inspect it for damage or deterioration. If there is any visible deterioration (cracking, tears, etc), do not use your mask and seek replacement part(s) from your healthcare provider.
 - ii. Inspect the exhaust holes on the Elbow **E**. Ensure air is flowing through these holes and do not use if blocked.
 - iii. Make sure path of air in mask is free from objects that may prevent breathing.

NOTE: Failure to follow the operating instructions above may compromise the performance and safety of the mask.

i MASK PARTS


Refer to the mask diagram on the inside cover of this guide.

- | | |
|-------------------------|--------------------|
| A Silicone Seal | E Elbow |
| B Mask Frame | F Swivel |
| C Headgear | G Diffuser |
| D Headgear Clips | H Mask Tube |

Note: The Elbow **E** and Mask Tube **H** are not intended to be removed from the mask by the home user. If removed unintentionally, refer to technical information page for assembly.


ii FITTING YOUR MASK

1. Remove all packaging before using your mask.
2. Hold the Silicone Seal **A** away from your nose, use your fingers to carefully spread open the soft side wings of your Silicone Seal **A** and then guide the two openings into your nostrils.
3. Hold the Mask Frame **B** in position with one hand and pull the Headgear **C** over your head with the other hand, positioning the straps above the ears as per the fitting images.
4. Gently tighten or loosen the Headgear straps on the sides of the Headgear **C** until you achieve a comfortable seal.

Useful fitting tips:

- If you are experiencing any leaks after fitting your mask, you may need to lift the Silicone Seal **A** off your nose and reposition it again so that the two openings on the Silicone Seal **A** fit securely into your nostrils.
- Ensure the mask headgear straps are not fastened too tight.
- Ensure the two silicone openings are not inverted into the Seal **A** and can be placed into each nostril.

iii DISASSEMBLY FOR CLEANING

1. Unhook the Headgear Clips **D** from the Mask Frame **B**.
2. Remove the Silicone Seal **A** from the Mask Frame **B**.
3. Unclip the Diffuser **G** from the Elbow **E**.
4. Remove the Swivel **F** from the Mask Tube **H**.

iv CLEANING YOUR MASK AT HOME

Wash your mask (excluding the Headgear **C** and Diffuser **G**) after each use:

1. Hand-wash your mask in soap dissolved in lukewarm water. Do not soak for more than 10 minutes.
2. Rinse thoroughly with fresh water, ensuring that all soap residue has been removed.
3. Inspect all parts to ensure they are visually clean. If necessary, repeat the cleaning steps until all parts are visually clean.
4. Leave all parts to dry out of direct sunlight before reassembling.

Wash your Headgear **C** and Diffuser cover **G** after 7 days of use:

1. Hand-wash your Headgear **C** in soap dissolved in lukewarm water. Do not soak your Headgear for more than 10 minutes.
2. Gently hand wash the Diffuser **G** in soap dissolved in lukewarm water, avoid scrubbing the material. Do not soak for more than 10 minutes.
3. Rinse thoroughly with fresh water, ensuring that all soap residue has been removed.
4. Inspect all parts to ensure they are visually clean. If necessary, repeat the cleaning steps until all parts are visually clean.
5. Leave all parts to dry out of direct sunlight before reassembling.

Cleaning Cautions:

- Do not clean your mask with products containing alcohol, anti-bacterial agents, antiseptic, bleach, chlorine or moisturizer.
- Do not clean your mask in a dishwasher.
- Do not store your mask in direct sunlight.

The above actions may deteriorate or damage the mask and shorten its life. If there is any visible deterioration do not use your mask and seek replacement part(s).


V MASK ASSEMBLY

1. Push the Silicone Seal **A** onto the Mask Frame **B** in the orientation shown in the below image until a confirmation click is heard.
2. Attach the Headgear Clips **D** on the Headgear **C** to the Mask Frame **B**.
3. Attach the Diffuser **G** to the Elbow **E**.
4. Push the Swivel **F** onto the Mask tube.


Mask Assembly Tip:

- Check the Swivel **F** is not stuck in the CPAP tube.


1


2


3

4

SYMBOL DEFINITIONS

	Catalogue number
	Manufacturer
	Authorized European Representative
	Batch code
	Temperature limitation
	Date of manufacture
	Consult instructions for use
	Caution
	Conforms with medical device directive 93/42/EEC
Rx only	Prescription only
	Extra small to small seal
	Medium to large seal

OXYGEN/PRESSURE PORT CONNECTOR

If pressure readings and/or additional oxygen are required, an Oxygen/Pressure Port Connector is available (REF 900HC452). Please contact your healthcare provider.


TECHNICAL SPECIFICATIONS

F&P BREVIDA EXHAUST FLOW

The F&P Brevida Nasal Pillows Mask has a unique Diffuser that is designed to diffuse air that is expelled from the mask. This reduces the draft. The Diffuser does not significantly change the amount of air flowing from the mask, therefore the mask can be used with or without the Diffuser.

Pressure (cmH ₂ O)	4	5	7	9	11	13	15	17	19	21	23	25
Flow with Diffuser (L/min)	21	24	29	33	37	41	44	47	50	53	56	58
Flow without Diffuser (L/min)	25	27	32	37	41	44	48	51	54	56	59	62

Due to manufacturing variations, the exhaust flow rates may vary from the nominal values shown above.

RESISTANCE TO FLOW

Pressure drop through XS-S mask with Diffuser @ 50 L/min: 1.7cmH₂O ± 10%
 Pressure drop through M-L mask with Diffuser @ 50 L/min: 1.0cmH₂O ± 10%
 Pressure drop through XS-S mask with Diffuser @ 100 L/min: 6.3cmH₂O ± 11%
 Pressure drop through M-L mask with Diffuser @ 100 L/min: 3.8cmH₂O ± 11%

MASK DEAD SPACE

XS-S: 29 cc

M-L: 33 cc

SOUND

The sound power level of mask with Diffuser is 25.4 dBA, with uncertainty 2.5 dBA.
 The sound pressure level of the mask with Diffuser is 17.5 dBA, with uncertainty 2.5 dBA.
 The sound power level of mask without Diffuser is 30.6 dBA, with uncertainty 2.5 dBA.
 The sound pressure level of the mask without Diffuser is 22.6 dBA, with uncertainty 2.5 dBA.

The F&P Brevida Nasal Pillows Mask and accessories are not made with natural rubber latex.

For single patient reuse, the mask should be replaced no later than 12 months after the date of first use. If any visible deterioration of a mask system component is apparent (for example cracking, discoloration, or seal tears), the component should be discarded and replaced.

MULTI-PATIENT CLEANING INSTRUCTIONS

For the USA only, BREISA, BREIMA and BREISMA are single patient use only, therefore multi-patient cleaning instructions do not apply.

Your F&P Brevida Nasal Pillows Mask is provided patient-ready.

To reuse between patients please follow the guidelines below or see the Fisher & Paykel Healthcare website: www.fphcare.com/disinfection

If you do not have internet access, please contact your Fisher & Paykel Healthcare representative.

HEADGEAR (C AND D)

1. Hand-wash the Headgear (C and D) in soap dissolved in lukewarm water. Air-dry out of direct sunlight.

DIFFUSER (G)

2. The Diffuser (G) is for single patient use. Replace between patients.

ALL OTHER PARTS

3. Soak all reprocessable parts of the mask in a solution of mildly alkali, anionic detergent e.g. Alconox according to the manufacturer's instructions (at 1% concentration and 50°C for 10mins), ensuring that no air bubbles are present on the device surface. Using a soft non-metallic brush (e.g. medium-hard toothbrush*), vigorously brush the mask parts until visibly clean but not for less than 10 seconds. Pay close attention to crevices and cannulations, such as the join between the silicone seal and plastic housing and the socket between the elbow and frame, and elbow and tube. Rinse by submerging in 5 litres of demineralised water and agitating for at least 10 seconds. Repeat rinse.
*An additional thin bottle brush is required to access the inside of the Brevida Nasal Pillows Mask tube.

4. HIGH-LEVEL THERMAL DISINFECTION	HIGH-LEVEL CHEMICAL DISINFECTION	STERILIZATION
Pasteurize at 75 °C for 30 mins, 80 °C for 10 mins or 90 °C for 1 min.	CIDEX OPA/METRICIDE OPA PLUS/RAPICIDE OPA/28: Immerse at 20 °C (68 °F) for 12 minutes then rinse thoroughly according to the manufacturer's rinsing instructions. (Replace seal. Replace tube if discolours)	Utilize a STERRAD NX, 100S or 100NX Sterilization System according to the manufacturer's instructions, using the pouch method.

Before each use, inspect the mask for deterioration. If any visible deterioration of a mask component is apparent, the mask component discarded and replaced. Some discoloration may occur after reprocessing and is acceptable.

To ensure traceability, mask components belonging to an individual mask must be reassembled back into the same mask after cleaning.

The F&P Brevida Nasal Pillows Mask has been validated to withstand 20 cycles of the above parameters. Failure to properly clean this device may result in inadequate disinfection. Fisher & Paykel Healthcare does not recommend any deviations from the recommended method of reprocessing.

ADVANCED DISASSEMBLY/ASSEMBLY

Please refer to the diagrams at the back of this guide.

- 1 **To disassemble Mask Tube (H) from Elbow (E):**
Pull the rigid connection on the Mask Tube (H) down from the bottom of the Elbow (E) until it disconnects.
- 2 **To assemble Mask Tube (H) onto Elbow (E):**
Push the rigid connection on the Mask Tube (H) onto from the bottom of the Elbow (E) until it clicks into position.
- 3 **To disassemble Elbow (E) from Mask Frame (B):**
 - Hold the Mask Frame (B) securely with one hand and position the Elbow (E) in its upright orientation.
 - Pull the bottom of the Elbow (E) upward from the Mask Frame (B) until it dislodges completely.
- 4 **To assemble Elbow (E) onto Mask Frame (B):**
 - Hold the Mask Frame (B) securely with one hand and position the top surface of the Elbow (E) ball inside the top edge of the Mask Frame (B) socket.
 - Apply downward pressure on the Elbow (E) until the bottom surface of the Elbow (E) enters the socket.
 - Ensure the Elbow (E) can rotate freely once inserted.

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